




Animal Health Product

Giving Set & infusion Consumables

2021



Primary Set..... 2~5
Extension Set.....6
Accessories.....7



As a manufacturer, BQ+'s choice of innovative infusion consumables including pump compatible infusion set, extension set and accessories, etc.

We would help deliver comprehensive solution for all your fluid therapy needs. Welcome to BQ+ intravenous fluid therapy online catalogue. We are well aware that fluid therapy forms a vital component of veterinary medicine. The role of a properly managed fluid therapy plan is increasingly recognised as a major contributing factor to the successful outcome in the treatment of many medical conditions.

Phone:(86)21-5774 3953
www.bqplusmedical.com

Intravenous Fluid Therapy

Primary set Easy-Order Table

| Product | Type | Micro | Macro | Big Bore | Light Protection | Unit/box | BQ+ Code | |
|--|--|-------|-------|----------|------------------|----------|--|-------------------------------|
|  Series A | -Compatible with Baxter Pump | | | | | |  Series B | -Compatible with Hospira Pump |
| Application Pump Giving Sets | | | | | | | Series A | Series B |
| Pump Giving Set | With 2 injection sites Full pump tubing | | ✓ | | ✓ | 50 | A001 | B001 |
| | With 2 needle free valves Full pump tubing | | ✓ | | ✓ | | A002 | B002 |
| | With 2 injection sites Partly pump tubing | | ✓ | | ✓ | | A003 | B003 |
| Measured Volume Gravity Infusion Sets | | | | | | | | |
| Dial Flow Type | With Precision Dial Flow 5-250ml adjustable | | ✓ | | ✓ | 40 | GD001 | |
| Burette Type | With Burette from 0-150ml Auto shut off valve | | ✓ | | ✓ | 100 | GB001 | |


Extension Set Easy-Order Table

| Product | Type | Micro | Macro | Big Bore | Light Protection | Unit/box | BQ Series # | |
|----------------------|-----------------------------|-------|-------|----------|------------------|----------|--------------|--------------|
| Extension Set | | | | | | | Micro | Macro |
| Dial Flow Type | With In line flow regulator | | ✓ | | | 50 | ED001 | |
| Coiled Set | With releaseable tie | | ✓ | | | | EC001 | |
| J Loop Set | With anti kink Design | | ✓ | | | | EJ001 | |
| T injection Site Set | With/o Needle free valve | ✓ | ✓ | | | | ET001 | ET002 |

Accessories Easy-Order Table

| Product | Type | Micro | Macro | Big Bore | Light Protection | Unit/box | BQ Series # | |
|--------------------------|--------------------------|-------|-------|----------|------------------|----------|-------------|--|
| Accessory | | | | | | | | |
| Transfer set | With double spike | | | ✓ | | 50 | T001 | |
| Needle free Valve | Male to Female luer ends | | | | ✓ | 100 | NFV-L16 | |
| I.V. Bag Access Device | with needle free valve | | | | | 100 | ST18-NF-01 | |
| Needle free Vial Adapter | with 0.22um filter | | | | | 100 | UVA-NF-01 | |

Primary set–Series A–Compatible with Baxter Pump




| BQ+ Code | Product Description | Size | Units/case |
|----------|---------------------|------|------------|
|----------|---------------------|------|------------|

| | | | |
|------|---|---------|----|
| A001 | <p>Pump Giving Set</p> <p>Vented Drip Chamber, 10 drops/mL, full pump tubing Series A slide clamp, roller clamp, 2 injection ports, Rotating LL.</p> | 78 inch | 50 |
|------|---|---------|----|




| | | | |
|------|--|---------|----|
| A002 | <p>Pump Giving Set</p> <p>Vented Drip Chamber, 10 drops/mL, full pump tubing Series A slide clamp, roller clamp, 2 needle free valves, Rotating LL.</p> | 78 inch | 50 |
|------|--|---------|----|



| | | | |
|------|--|---------|----|
| A003 | <p>Pump Giving Set</p> <p>Vented Drip Chamber, 10 drops/mL, part pump tubing Series A slide clamp, roller clamp, 2 needle free valves, Rotating LL.</p> | 78 inch | 50 |
|------|--|---------|----|



Primary set–Series B–Compatible with Hospira Pump



| BQ+ Code | Product Description | Size | Units/case |
|----------|---------------------|------|------------|
|----------|---------------------|------|------------|

| | | | |
|------|---|---------|----|
| B001 | Pump Giving Set Vented Drip Chamber, 10 drops/mL, full pump tubing Series B slide clamp, roller clamp, 2 injection ports, Rotating LL. | 78 inch | 50 |
|------|---|---------|----|




| | | | |
|------|--|---------|----|
| B002 | Pump Giving Set Vented Drip Chamber, 10 drops/mL, full pump tubing Series B slide clamp, roller clamp, 2 needle free valves, Rotating LL. | 78 inch | 50 |
|------|--|---------|----|



| | | | |
|------|--|---------|----|
| B003 | Pump Giving Set Vented Drip Chamber, 10 drops/mL, part pump tubing Series B slide clamp, roller clamp, 2 needle free valves, Rotating LL. | 78 inch | 50 |
|------|--|---------|----|



Primary set–Measured Volume Gravity Infusion Sets



| BQ+ Code | Product Description | Size | Units/case |
|----------|---------------------|------|------------|
|----------|---------------------|------|------------|

GD001

Dial Flow Type

Vented Drip Chamber, 20 drops/mL,
Precision flow regulator 5-300ml, Rotating LL.

78 inch

50




GB001

Burette



Vented Drip Chamber, 20 drops/mL,
With auto shut off valve, from 0-150ml, Rotating LL.
Series B slide clamp, roller clamp, 2 needle free valves, Rotating LL.






78 inch

50





Extension Set



| BQ+ Code | Product Description | Size | Units/case |
|----------|--|---------|------------|
| ED001 | <p>Dial Flow Type</p> <p>Winged connector, in-line flow regulator, injection port, lockable slip luer.</p>  | 78 inch | 50 |
| EC002 | <p>Coiled Set</p> <p>Extra Thickness Tubing, Male & Female connector With releasable tie to fix on collar</p>  | 78 inch | 50 |
| EJ001 | <p>J Loop Set</p> <p>Extra Thickness Tubing, Male & Female connector With J Loop design to fix tubing</p>  | 78 inch | 50 |
| ET001 | <p>T injection site set</p> <p>Micro Bore tubing, T injection site luer lock</p>  | 78 inch | 50 |
| ET002 | <p>T injection site set</p> <p>Micro Bore tubing, T injection site luer lock</p>  | 78 inch | 50 |

Accessories



| BQ+ Code | Product Description | Size | Units/case |
|----------|---------------------|------|------------|
|----------|---------------------|------|------------|


| | | | |
|------|--|---------|----|
| T001 | Transfer Set Big bore transfer set for big animal | 78 inch | 50 |
|------|--|---------|----|



| | | | |
|---------|--|---|----|
| NFV-L16 | Needle Free Valve Male to Female luer ends, single packed | / | 50 |
|---------|--|---|----|



| | | | |
|------------|---|---|----|
| ST18-NF-01 | I.V. Bag Access Device Needle-free injection of additives into fluid bag | / | 50 |
|------------|---|---|----|



| | | | |
|-----------|---|---|----|
| UVA-NF-01 | Needle free Vial Adapter For needle-free aspiration or | / | 50 |
|-----------|---|---|----|

